



Denville XL3000i pipettes



- Classic profile
- Continuously adjustable
- Adjustable tip ejector is easily removed
- Comfortable, ergonomic hand grip
- Liquid end is autoclavable
- Seven different models covering the range from 0.1µl to 5000µl

loose and racked tips are available for Denville pipettes
 please call for details

DENVILLE 
SCIENTIFIC INC.



also available for 3x pipettes

Seven models of Denville XL3000i™ Pipettes cover all volumes from 0.1µl up to 5000µl. A continuously adjustable volumeter with digital readout allows simple and accurate dispensing. The robust construction along with the low thermal coefficient of the body of each XL3000i™ will prevent hand heat affecting sample measurements and reproducibility, even in prolonged usage. Being lighter than many other similar pipettes, they are even more comfortable to hold and operate for extended periods - even in the smallest hand. This feature will also reduce the incidence of operator fatigue and Repetitive Strain Injury.

The XL3000i™ requires minimal maintenance. Its precise, self-locking stainless steel micrometer accurately adjusts the stroke of the polished, acid-resistant stainless steel piston. Each piston is machined to very fine tolerances which, in conjunction with its self-lubricating O-ring and PTFE seal, minimises the need for attention. Using a series of height adapters, the tip ejector position on the shaft is adjustable to allow the pipette's use with virtually all brands of tips. The user may also choose from either of the two supplied ejector button styles to suit their preferred handling position.

Each model is supplied with its own individual certificate of calibration, as a guarantee of the unit's quality.

try before you buy†

1 EJECTOR BUTTON
two styles to maximise individual comfort

2 EASY-TURN
volume adjustment wheel

2 DIGITAL
volume display

SLIM HANDLE
for optimum comfort

LIGHT WEIGHT
- up to 20% lighter than many similar models

PLASTIC GRIP
for easier tip ejector removal

3 AIRTIGHT SEAL
PTFE / O-ring

4 ADJUSTABLE EJECTOR
to allow a wider choice of tips

1 Conventional holding style - standard ejector button

1 Side holding style - large ejector button for increased comfort

2 Easy read digital display

3 Comprehensive spares

4 Easy removal of tip ejector

† try a **Denville pipette** for **7 days**. If not completely satisfied, return it in the packaging provided. (UK only)

Ordering Information

Description	Accuracy %	Coeff. Variation %	Cat. No.
0.1 to 2µl	±12 to 1.5	±6.0 to 0.7	P-3950-2
0.5 to 10µl	±3.6 to 0.5	±3.6 to 0.3	P-3950-10
2 to 20µl	±4.0 to 0.8	±3.0 to 0.4	P-3950-20
10 to 100µl	±1.6 to 0.8	better than ±1%	P-3950-100
20 to 200µl	±1.2 to 0.6	better than ±1%	P-3950-200
100 to 1000µl	±0.9 to 0.6	better than ±1%	P-3950-1000
1 to 5ml	±0.6 to 0.5	better than ±1%	P-3950-5000
1 to 10ml	±2.5 to 0.5	better than ±1%	P-3950-10ML
Pipette Stand, 4-position	BC-4020	Pipette Stand, 3-position	BC-4021